SKIN MARKERS

REF	PS3150, PS3150L, PS3151, PS3158, PS3150NR
	en – Skin Marker es – Rotulador cutáneo pt – Marcador de Pele de – Hautmarkerierungsstift fr – Crayon dermographique it – Pennarello dermografico cs – Kožní značkováč da – Tusch til hud sv – Hudmarkeringspenna no – Hudmarkør et – Nahamarker hu – Bőrfelölő lt – Odos žymeklis pl – Marker do skóry tr – Deri İşaretleyici ro – Marker de piele

Purple Surgical International Limited/Purple Surgical UK Limited
 2 Chestnut House, Farm Close, Shenley, Herts WD7 9AD
 England T: +44(0)1923 839333 | F: +44(0)1923 839444
 E: mail@purplesurgical.com | W: www.purplesurgical.com

Advena Limited
 Tower Business Centre
 2nd Floor Tower Street
 Swatar, BKR 4013
 Malta

Medenvoy Switzerland
 Gotthardstrasse 28
 6302 Zug
 Switzerland

Medenvoy Global BV
 Prinses Margrietplantsoen
 33 – Suite 123
 2595 AM The Hague
 The Netherlands

Medenvoy Switzerland
 Gotthardstrasse 28
 6302 Zug
 Switzerland

Skin Marker

The following should be thoroughly read before use.

IMPORTANT NOTES

This document is designed to provide instructions on how this product/s is used. It is not a reference to clinical techniques and does not replace local or national guidance or protocols.

INTENDED USER

This product should only be used by suitably trained healthcare professionals within a healthcare establishment.

INTENDED PURPOSE

A surgical pen designed specifically for clinical markings on the skin to promote correct site surgery.

INDICATIONS FOR USE

Use is governed by the clinical markings on the skin to promote correct site surgery.

Functional uses are:

- Ensure surgery is performed on the correct site.
- Mark an important internal anatomical feature identified during diagnostic imaging.
- Target non-invasive treatments.
- Plan the surgical approach.
- Ensure natural re-alignment of skin flaps during wound closure.

INTENDED PATIENT POPULATION

There are no patient population limitations for these devices.

PRODUCT CHARACTERISTICS

The Purple Surgical Skin Marker allows for the precise marking on skin prior to surgery.

HOW SUPPLIED

Supplied sterile, ready for single patient use.

SHELF LIFE

4 years from date of sterilisation - Use by date as denoted on labelling in format YYYY-MM

OTHER DEVICES REQUIRED FOR USE AND INTEROPERABILITY REQUIREMENTS

The Purple Surgical Skin Marker is not designed to connect to other devices and has no interoperability requirements.

CONTRAINDICATIONS

Do not use on broken / breached skin. Only use on intact skin.

PRECAUTIONS / WARNINGS / RESIDUAL RISKS

- (1) DO NOT use if package is received opened or damaged.
- (2) Procedures should be performed only by clinicians with adequate training and knowledge. Medical literature should be consulted for techniques, hazards, contraindications and complications prior to the procedure.

POSSIBLE ADVERSE REACTIONS

No known issues.

SINGLE USE PRECAUTION

These device(s) are designed and sold for single patient use only.

Re-processing and/or re-sterilisation is not permitted. The effects of unauthorised reprocessing or re-sterilisation can result in the following complications:

1. Infection due to cross contamination as a result of ineffective re-processing/re-sterilisation.
2. Mechanical fatigue and performance failure, due to the effects of the re-processing / re-sterilisation method.

DIRECTIONS FOR USE

- (1) Visually inspect the sterile barrier system pouch prior to its opening so as to ensure it is not already open or compromised. If a sterile barrier failure is found, the contained product shall be considered non-sterile and it shall not be used.
- (2) Using sterile technique, remove the product/s from its packaging. **NOTE:** Begin opening the sterile pouch at the pull tab seal end.
- (3) Apply the single-use Skin Marker to the area of dry skin intended for surgery. Allow to dry for 60 seconds before prepping the skin with iodine.

DISPOSAL

Discard after single patient use, DO NOT re-sterilise.

This product is to be disposed of as controlled medical waste according to national guidelines.

SERIOUS INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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