




PRIME MEDICAL llc

CLINICALLY DRIVEN INNOVATION



SURGICAL PATIENT POSITIONERS AND ACCESSORIES



CLINICALLY DRIVEN INNOVATION

Prime Medical, LLC was formed in 2012 by a group of individuals with diverse backgrounds in medical and surgical devices who have a passion for innovation with patient focus. Prime Medical's products are developed through collaboration with a variety of clinicians, including surgeons, anesthesiologists, and O.R. nurses. Many of the devices in Prime Medical's arsenal solve patient positioning and safety problems not currently being addressed in the Operating Room. We welcome ideas from all surgical services and applications, as we are constantly developing clinically driven innovation in a wide variety of clinical areas.

Prime Medical, LLC is based in Knoxville, TN and each product in this catalog is proudly made in the U.S.A.

SURGICAL PATIENT POSITIONERS AND ACCESSORIES

Pages 3-10

PrimePad Anti-Skid O.R. Pad Series

Page 11-15

Anesthesia Positioners

Pages 16-20

Orthopaedic and Spine Positioning Accessories

Pages 21-23

Pediatric Positioners

Pages 24-27

Patient Safety Accessories



PRIME MEDICAL SURGICAL SOLUTIONS

PrimePad Anti-Skid O.R. Pad Series

- For use in any procedure that requires Trendelenburg, Reverse Trendelenburg, or lateral positioning.
- Provides superior support and stability with a next-generation, innovative, high-performance polyurethane foam pad.
- PrimePad material is made with microscopic air capsules that offer immediate rebound, anti-skid properties, enhanced breathability and increased patient support.
- PrimePad material is not temperature sensitive and dissipates heat faster than memory foam.
- PrimePad material provides up to 4 times greater peak pressure reduction compared with memory foam.
- Single patient use due to the dissipation of static properties with each use.
- Single patient use to aid in infection prevention.
- Compatible with patient warming devices.
- Single-packed in compressed format which lowers shipping costs and reduces storage space requirements.
- Available in kits that include: Patient Body Straps, Disposable Draw Sheets, Arm Tucking/Protection Sleeves, and Head Cradles.
- Available in Wing format with integrated arm tucking and protection.

PrimePad-1



Case Pack: 10

PrimePad-1S



Case Pack: 10



Instructions for Use

- First, ensure the O.R. Mattress is secured to the table.
- O.R. Mattress must be clean and dry. Place pad on the table, aligning the pad and mattress perineal cutouts. PrimePads are designed to provide maximum O.R. Mattress coverage.
- Nothing should be placed between O.R. mattress and PrimePad.
- MAXIMUM patient contact with the pad is key to clinical efficacy.
- Ensure the Sacrum, Torso and Shoulders are in DIRECT contact with the pad.
- Minimize lift/draw sheet width to maximize patient contact with the pad.
- Prefer foam or gel donuts over pillows for intra-operative head positioning, as patients' shoulders often rest on pillow instead of the pad. If using pillows, ensure shoulders are in contact with pad.
- Pads with Table Rail Straps Only: The built-in table rail straps should be on the bottom side of the pad. Slide the straps downward between the rail and the table, bring them back up and around the outside of the rail, and affix the strap back upon itself securely. Discard after use.

Not made with natural rubber latex



PrimePad Anti-Skid O.R. Pad Series

PrimePad-1: Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 2.54 cm.

PrimePad-1S: Anti-Skid Pad w/perineal cutout w/rail straps, 112 cm x 58 cm x 2.54 cm.

PrimePad-2: 5.08 cm Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 5.08 cm.

PrimePad-2S: 5.08 cm Anti-Skid Pad w/perineal cutout w/rail straps, 112 cm x 58 cm x 5.08 cm.

PrimePad-3: Ramped 2.54 cm to 5.08 cm Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 5.08 cm.

PrimePad-3S: Ramped 2.54 cm to 5.08 cm Anti-Skid Pad w/perineal cutout, rail straps, 112 cm x 58 cm x 5.08 cm.

PrimePad-W1: Anti-Skid Pad w/integrated 61 cm x 61 cm arm tucking/protection, perineal cutout, rail straps, 112 cm x 58 cm x 2.54 cm.

PrimePad-W2: Anti-Skid Pad w/integrated 61 cm x 41 cm arm tucking/protection, perineal cutout, rail straps, 112 cm x 58 cm x 2.54 cm.

PrimePad-K1: Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, 112 cm x 58 cm x 2.54 cm.

PrimePad-K2: Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, disposable draw sheet, 112 cm x 58 cm x 2.54 cm.

PrimePad-K3: Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, arm tucking sleeve pair, 112 cm x 58 cm x 2.54 cm.

PrimePad-K4: Anti-Skid Pad w/perineal cutout, head/shoulder positioner, single body strap set, arm board pair, 112 cm x 58 cm x 2.54 cm.

PrimePad-K5: 5.08 cm Anti-Skid Pad w/perineal cutout, head/shoulder positioner, single body strap set, arm tucking sleeve pair, 112 cm x 58 cm x 5.08 cm.

PrimePad-K6: Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, disposable draw sheet, arm tucking sleeve pair, 112 cm x 58 cm x 2.54 cm.

PrimePad-K7: Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, disposable draw sheet, arm tucking sleeve pair, head cradle, 112 cm x 58 cm x 2.54 cm.

PrimePad-K8: Anti-Skid Pad w/integrated 61 cm x 41 cm arm tucking/protection, perineal cutout, rail straps, single body strap set, disposable draw sheet, head cradle, 112 cm x 58 cm x 2.54 cm.

PrimePad-K9: Anti-Skid Pad w/perineal cutout w/rail straps, disposable draw sheet, 112 cm x 58 cm x 2.54 cm.

PrimePad-K10: 5.08 cm Anti-Skid Pad w/perineal cutout w/rail straps, single body strap set, disposable draw sheet, 112 cm x 58 cm x 5.08 cm.

PrimePad-PED: Pediatric Anti-Skid Pad, 61 cm x 51 cm x 2.54 cm.

PrimePad-SQ1: Square positioning foam, 30 cm x 30 cm x 2.54 cm.

PrimePad-SQ2: 5.08 cm Square positioning foam, 30 cm x 30 cm x 5.08 cm.

PrimePad-BBO: Bean Bag Overlay, 102 cm x 97 cm x 2.54 cm.

PrimePad-CVP: Cardiovascular Pad, 183 cm x 58 cm x 2.54 cm.



PRIME MEDICAL SURGICAL SOLUTIONS

PrimePad Anti-Skid O.R. Pad with Integrated Arm Tucking and Protection

-  Provides all of the features and benefits of our PrimePad Series, with the added feature of integrated arm tucking and protection Wings.
-  The integrated arm tucking and protection offers patient safety and clinical staff ease of use and convenience in patient positioning for lithotomy/steep Trendelenburg.
-  The Wing material can be used to either fully encapsulate the patients' arms, or they can be used to simply pad the patients' arms, depending upon hospital protocol.
-  Any excess or unnecessary Wing material can be removed by simply tearing it, and it easily tears in a linear plane.
-  Allows the Anesthesia provider easy access to the IV site, and Wing material allows for channeling of IV tubing and pulse oximetry cables by tearing a window or slit in the Wing material.
-  Eliminates the need to use a draw sheet for arm tucking and allows for 100% patient contact with the pad.
-  2 options with different widths in arm tucking/ protection material:

PrimePad-W1 Wing dimensions:
61 cm long x 61 cm wide

PrimePad-W2 Wing dimensions:
61 cm long x 41 cm wide

PrimePad-W2



Case Pack: 4



Providing a Better Way

The PrimePad is a new innovation in anti-skid patient positioning, offering a next-generation, innovative, high-performance polyurethane foam pad. PrimePad offers excellent anti-skid properties, pressure reduction and patient support, providing confidence to the surgeon and staff that the patient will remain immobile.

Instructions for Use

Follow the Instructions For Use for the PrimePad Series of Anti-Skid Pads for installation of the pad on the O.R. Mattress.

Once the patient is properly positioned on the PrimePad, a draw sheet is not required to tuck and protect the arms. To allow for 100% patient-to-pad contact, remove the draw sheet and utilize the Wing material to tuck/protect the arms.

To utilize the Wing material for arm tucking and protection, bring the full length of the pad to the patient's side to determine how much length is needed. If excess Wing material needs to be removed, simply tear off the un-needed material. The material tears in linear plane, and tearing reduces the need for scissors.

One method for utilizing the Wing material is to have a team member on each side of patient. One team member will slightly "log roll" the patient toward them, while the other team member wraps the patient's arm in the Wing and then places the end of the Wing under the patient's torso. Repeat the same steps on the opposite side.

The PrimePad material's flexibility and high coefficient of static friction provides secure arm tucking and protection without the need for external accessories such as Velcro, straps and draw/lift sheets.

Not made with natural rubber latex



For PrimePads with Integrated Arm Tucking and Protection:





PrimePad Series

All PrimePads are folded, compressed and vacuum sealed for compact storage and shipping efficiency.



PrimePad-1
Case Pack: 10



Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 2.54 cm

PrimePad-1S
Case Pack: 10



Anti-Skid Pad w/perineal cutout w/rail straps, 112 cm x 58cm x 2.54 cm

PrimePad-2
Case Pack: 4



5.08 cm Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 5.08 cm

PrimePad-2S
Case Pack: 4



5.08 cm Anti-Skid Pad w/perineal cutout w/rail straps, 112 cm x 58 cm x 5.08 cm

PrimePad-3
Case Pack: 4



Ramped 2.54 cm to 5.08 cm Anti-Skid Pad w/perineal cutout, 112 cm x 58 cm x 5.08 cm

PrimePad-3S
Case Pack: 4



Ramped 2.54 cm to 5.08 cm Anti-Skid Pad w/perineal cutout, rail straps, 112 cm x 58 cm x 5.08cm

PrimePad Series

PrimePad-W1

Case Pack: 4



Anti-Skid Pad w/integrated
61 cm x 61 cm arm tucking/
protection, perineal
cutout, rail straps, 112 cm x
58 cm x 2.54 cm

PrimePad-W2

Case Pack: 4



Anti-Skid Pad w/integrated
61 cm x 41 cm arm tucking/
protection, perineal
cutout, rail straps, 112 cm x
58 cm x 2.54 cm

PrimePad-K1


Case Pack: 10



Anti-Skid Pad w/perineal
cutout w/rail straps,
single body strap set,
112 cm x 58 cm x
2.54 cm

PrimePad-K2

Case Pack: 10



Anti-Skid Pad w/perineal
cutout w/rail straps, single
body strap set, disposable
draw sheet, 112 cm x 58
cm x 2.54 cm

PrimePad-K3

Case Pack: 3



Anti-Skid Pad w/perineal
cutout w/rail straps, single
body strap set, arm tucking
sleeve pair, 112 cm x 58 cm
x 2.54 cm

PrimePad-K4

Case Pack: 3



Anti-Skid Pad w/perineal
cutout, head/shoulder
positioner, single body
strap set, arm board
pair, 112 cm x 58 cm x
2.54 cm

PrimePad Series

PrimePad-K5

Case Pack: 3



5.08 cm Anti-Skid Pad w/
perineal cutout, head/
shoulder positioner, set of
two body straps, arm
tucking sleeve pair, 112 cm
x 58 cm x 5.08 cm

PrimePad-K6

Case Pack: 4



Anti-Skid Pad w/perineal
cutout w/rail straps, single
body strap set, disposable
draw sheet, arm tucking
sleeve pair, 112 cm x 58
cm x 2.54 cm

PrimePad-K7

Case Pack: 3



Anti-Skid Pad w/perineal
cutout w/rail straps, single
body strap set, disposable
draw sheet, arm tucking
sleeve pair, head cradle,
112 cm x 58 cm x 2.54 cm

PrimePad-K8

Case Pack: 3



Anti-Skid Pad w/integrated 61
cm x 41 cm arm tucking/
protection, perineal cutout,
rail straps, single body strap
set, disposable draw sheet,
head cradle, 112 cm x 58 cm x
2.54 cm

PrimePad-K9


Case Pack: 10



Anti-Skid Pad w/perineal
cutout w/rail straps,
disposable draw sheet,
112 cm x 58 cm x
2.54 cm

PrimePad-K10

Case Pack: 4



5.08 cm Anti-Skid Pad w/
perineal cutout w/rail
straps, single body strap
set, disposable draw
sheet, 112 cm x 58 cm x
5.08 cm

PrimePad Series

PrimePad-PED

Case Pack: 10



Pediatric Anti-Skid
Pad, 61 cm x 51 cm x
2.54 cm

PrimePad-SQ1

Case Pack: 25



Square Positioning
Foam, 30 cm x 30 cm
x 2.54 cm

PrimePad-SQ2


Case Pack: 25



5.08 cm Square
Positioning Foam, 30
cm x 30 cm x 5.08 cm

PrimePad-BBO

Case Pack: 6



Bean Bag
Overlay, 102 cm x
97 cm x 2.54 cm

PrimePad-CVP

Case Pack: 4



Cardiovascular
Pad, 183 cm x 58
cm x 2.54 cm




PRIME MEDICAL SURGICAL SOLUTIONS

Lateral Positioning Pillow

- Properly aligns the patient's cervical spine.
- Channel protects the eye while allowing constant monitoring.
- Center cutout protects the ear.
- Pillow angle evenly distributes the weight of the head.
- Pillow height easily adjusts by removing one or both bottom layers to accommodate each patient.
- Provides consistency in lateral patient positioning.

LHP-100
US Patent # 8,516,639



Dimensions:

27 cm x 25 cm x 27 cm
(Height adjusts from 16.5 cm to 27 cm)

Case Pack: 8



Providing a Better Way

Consistency and standardization are the keys to patient safety when placing patients in the lateral position. Prime Medical's Lateral Positioning Pillow supports the patient's head for safe, reliable lateral positioning.

This two-part pillow has a firm, adjustable base that is anatomically angled to maintain the cervical spine in a neutral position. The soft-top layer conforms to the patient's face and evenly distributes the weight of the head throughout the lateral procedure.

The central cutout keeps the ear free of pressure. The patented eye channel allows the provider to constantly monitor the eye pressure to ensure it is free from pressure that could result in catastrophic injury.

The tear-away bottom layers afford the provider the ability to adjust the pillow to patients of varying sizes.

Not made with natural rubber latex



A] Improvised Positioning



B] Improper Alignment



C] Proper Positioning

Reliable Spine Alignment

Careful measurements of patients in the lateral position have demonstrated that the width of the shoulders compared to the hips and use of an axillary roll creates a natural angle of the thoracic and cervical spine relative to the horizontal plane of the surgical bed. Prime Medical's patented Lateral Positioning Pillow reproduces this angle consistently and reliably. Neutral thoracic and cervical spine alignment minimizes the risk of brachial plexus injuries.

Gone are the days of improvised support of the head in patients positioned laterally!




PRIME MEDICAL SURGICAL SOLUTIONS

Universal Supine/Prone Positioning Pillow

-  Properly aligns cervical spine in prone position.
-  Use for induction and positioning prone. Adjustable with easy tear-away layers.
-  Soft foam provides support and comfort.
-  Combines the supine induction head cradle and the prone pillow in a single positioner.
-  Trach slits on BOTH sides of the positioner.

UPSP100
US Patent # 9,138,087



Dimensions:
23 cm x 20 cm x 16.5 cm
(Height adjusts from 10 cm to 16.5 cm)

Case Pack: 10



Providing a Better Way

The bottom two layers can be removed so the pillow can be used in the supine position for induction. This all-in-one characteristic replaces the need for two pillows, and provides cost savings to both hospital and patient.

Instructions for Use


Key issues when positioning patients prone are maintaining a neutral cervical spine, and protecting the face and eyes. Prime Medical's Universal Prone Positioning Pillow addresses both of these. This two-part adjustable pillow allows the Anesthesia Provider to adjust the height of the pillow to accommodate patients of different sizes, as well as the use of varying chest support systems.




The ability to decrease the height of the pillow prevents hyper extension of the cervical spine when less bulky chest positioning systems are required. The firm base provides a stable, reliable platform while the soft inserts protect the face and eyes from pressure-related injuries. Not made with natural rubber latex



Now one pillow....




...does the job of two.




Cost Savings

The bottom two layers can be removed so the pillow can be used in the supine position for induction. This all-in-one characteristic replaces the need for two separate pillows, thus providing cost savings to both the hospital and patient.

Convenient
adjustment of
height with two
tear away layers.





CLINICALLY DRIVEN INNOVATION™


PRIME MEDICAL SURGICAL SOLUTIONS

Prone Head Positioner w/dual Trach Slit

The Prone Head Positioner (#PHP100) is indicated for all surgical procedures that require the Prone patient position.

- ▶ Measures 23 cm x 20 cm x 18 cm and offers 5.08 cm thick soft cosmetic foam top layer to greatly reduce the risk of skin injury, even in longer duration procedures.
- ▶ The Prone Head Positioner provides a 13 cm base layer with optimal density for patient support.
- ▶ Dual trach slits offer anesthesia providers maximum flexibility and ease of use.

PHP100



Case Pack : 10 eaches





PRIME MEDICAL SURGICAL SOLUTIONS

Universal Axillary/Positioning Roll

- ✚ Unique, adjustable single-patient-use foam positioner for a wide variety of surgical applications.
- ✚ Eliminates risk of cross-contamination and dermatological reaction to cleaning residue.
- ✚ Provides a superior alternative to rolled towels, bandages, blankets, sheets, I.V. bags, etc.
- ✚ Density of foam provides support with reduced risk of pressure-related tissue injury.
- ✚ Texture of foam decreases patient slippage.
- ✚ Flat bottom increases stability.
- ✚ Removable layers provide adjustability.
- ✚ Slot in the side will accommodate a 90-degree angle iron that can attach to a table rail clamp for added stability during knee procedures.



AAR100
US Patent # 9,161,642



Case Pack: 10

Dimensions:
32 cm long x 11 cm wide x 12 cm high
(2 cm tear-away layers)



Providing a Better Way

Consistency and standardization are the keys to patient safety when positioning patients for surgical procedures. Prime Medical's Universal Axillary/Positioning Roll safely and reliably supports the patient in multiple positions, bringing standardization to positioning for a wide variety of surgical procedures.

The design, density and texture of the foam provide stability and support, while the removable bottom layers allow the clinician to adapt the AAR100 to the patient and procedure. Use of the AAR100 eliminates the need for rolled towels, blankets, I.V. bags, etc., which are not positioning devices, or gel rolls which require cleaning between surgical cases.

Because the AAR100 is single-patient-use, the risk of cross-contamination from reusable positioners is eliminated. So, too, is the risk of a dermatological reaction to a cleaning solution residue. By removing the need for a cleaning protocol, patient safety is enhanced and OR time is saved.

The AAR100 is your consistent, safe and universal surgical positioning solution.




Horizontal Scapular Roll: Neck surgery - thyroidectomy, parathyroidectomy, radical neck dissection, internal jugular vein central line placement.


Not made with natural rubber latex



Application Examples



Axillary Roll: Lateral positioning - CV/thoracic, orthopaedic, nephrectomy




Vertical Scapular Roll: Median sternotomy; Subclavian central line placement



Foot/Ankle Prop: Podiatric surgery



Left Lateral Tilt: Cesarean section



PRIME MEDICAL SURGICAL SOLUTIONS


Chest Roll Pair

- Chest Rolls provide Standardization in Prone positioning requiring chest rolls by eliminating the need for rolled towels and other improvised positioning methods.
- Chest Rolls are adjustable in both length and height by removing the end and bottom portions to provide customized support for a wide variety of patient body habitus.
- Chest Rolls are single patient use to reduce the risk of cross-contamination and help increase room turnover by eliminating cleaning requirements.
- Chest Rolls offer ease of use for Clinical Staff by providing a lightweight, customizable and bacteriostatic positioning solution.

Not made with natural rubber latex



CR100



Case Pack: 7 Pair




PRIME MEDICAL SURGICAL SOLUTIONS

Bean Bag Overlay

- Provides skin protection and pressure reduction for Lateral positioning in conjunction with all types of bean bag devices.
- Made of the same next-generation, high-performance polyurethane foam as found in the PrimePad Series of Anti-Skid Pads for Trendelenburg, reverse-Trendelenburg and lateral positioning.
- Dimensions of 97 cm x 102 cm x 2.54 cm covers entire bean bag surface and allows maximum overlay coverage.
- Helps protect against nerve injury, compartment syndrome and pressure ulcers.
- Stored in compressed format which lowers shipping costs and reduces storage space requirements.

Not made with natural rubber latex

PrimePad - BBO



Case Pack: 6





PRIME MEDICAL SURGICAL SOLUTIONS

Arm Tucking Sleeve

-  Sleeves provide 360° padded protection for the arm and hands.
-  Sleeves allow easy access to IV sites for Anesthesia providers.
-  Sleeves easily accommodate blood pressure cuffs and IV tubing.
-  Sleeves are single-patient use and are latex free.

Not made with natural rubber latex



ATS100




Case Pack: 3 pair

Pediatric Positioners



Child Head Cradle

CHC100
Case Pack: 36



- ▣ Single patient use child head cradle for anesthesia positioning.
- ▣ Soft foam with appropriate density provides support and reduced pressure to the patient.

Not made with natural rubber latex



Pediatric Donut

PD100
Case Pack: 48



- ▣ Single patient use 18 cm donut for anesthesia positioning.
- ▣ Soft foam with appropriate density provides support and reduced pressure to the patient.

Not made with natural rubber latex





PRIME MEDICAL SURGICAL SOLUTIONS

SECURE360°

PEDIATRIC PRONE POSITIONER

- Provides 360-degree access and visualization of the infant's skull.
- Provides stable, reliable support throughout surgery.
- Significantly reduces patient positioning time.
- No need to reposition child during surgery.
- Intended for use with pediatric patients, specifically infants 6-18 months of age.



Case Pack: 1

PATENT US # 8,584,280

Providing a Better Way


Specifically designed for use during cranial reconstruction surgeries, the Secure360 Pediatric Prone Positioner reduces patient positioning time while providing stable support throughout surgery. Allowing 360 degree access and visualization of the infant's skull.

Secure360® Pediatric Prone Positioner is indicated for use during surgery for proper positioning of infants in the prone position while maintaining maximal visualization of the skull. It is designed to provide rapid, safe and stable prone positioning while specifically protecting the infant's eyes, neck, chest and extremities

Discard after use. Flame retardant. Silicone Free. Not made with natural rubber latex.

SECURE360°

PEDIATRIC PRONE POSITIONER



Heel Protector Pair

HP100

Case Pack: 36 pair




- Pair of single patient use convoluted foam heel protectors.
- Soft foam with appropriate density provides support and reduce pressure to the patient.
- Soft, wide foam straps for ease of use and patient safety.

Not made with natural rubber latex

Wheelchair Convoluted Foam Pad

WCF100

Case Pack: 12



- Single patient use foam pad for patient positioning in shoulder procedures utilizing a beach chair, as well as pressure reduction during wheelchair transportation.
- Soft foam with appropriate density provides support and reduced pressure to the patient.

- Dimensions: 9 cm x 41 cm x 46 cm.
- Not made with natural rubber latex

Jackson Table Rail Pads

JTRP100

Case Pack: 10 sets



- Single patient use set of rail pads (2-51 cm pads, 2-15 cm pads), compatible with most spine tables.
- Tear-off sections allow for adjustability in the length of the pads to accommodate various patient body habitus.

Not made with natural rubber latex

Adult Donut

AD9

Case Pack: 36



- Single patient use 23 cm donut for anesthesia positioning.
- Soft foam with appropriate density provides support and reduced pressure to the patient.

Not made with natural rubber latex

Rolled Foam Positioners

RFP100
Case Pack: 12




- Single patient use 43 cm x 12 cm foam roll for a variety of patient positioning options.
- Soft foam with appropriate density provides support and reduced pressure to the patient.

Not made with natural rubber latex

McGuire Positioner Pad System

MPPS100
Case Pack: 10 sets



- 4 piece single-patient use pad set for McGuire lateral positioning system.
- 2 anterior pads for adjustable arms, 1 posterior pad and 1 hip pad with peel-away adhesive for additional stability.

Not made with natural rubber latex

Laminectomy Arm Cradle Set

LAC100
Case Pack: 20 sets



- Set of 2 laminectomy arm cradles reduce pressure at bony prominences and decrease risk of pressure ulceration.
- Allows for optimal patient positioning.
- Also used in conjunction with arm-tucking during robotic-assisted Trendelenburg positioning cases for protection of ulnar nerve and stability of the arms.

Not made with natural rubber latex

Patient Safety/Body Straps

SS100 /SS50
Case Pack: 25 Sets/50 Sets




- Soft safety / body straps for use with O.R. table rails in a wide variety of patient positioning scenarios. (SS50 - single strap; SS100 - set of 2 straps).
- Single patient use straps accommodate all body types.

Not made with natural rubber latex

Eye Center Pillow

ECP 100
Case Pack: 10




Not made with
natural rubber
latex



- Single patient use head positioner for ophthalmological procedures.
- Provides stability for the patient's head and uses a soft forehead strap for additional security and comfort.
- Velcro material on bottom of pillow attaches to head rest on procedure table for optimal patient positioning and comfort.

Adult Head Cradle

HC100
Case Pack: 36



- Single patient use head cradle for anesthesia positioning.
 - Soft foam with appropriate density provides support and reduced pressure to the patient.
- Not made with natural rubber latex

Perineal Post


PPOST100
Case Pack: 40



- Relieves pressure in the perineal area during orthopedic fracture cases requiring mechanical traction.
 - Single patient use pad provides ease of use and reduced risk of cross contamination.
 - 30 cm x 9 cm
- Not made with natural rubber latex

Perineal Post

PPOST6
Case Pack: 12



- Relieves pressure in the perineal area during orthopedic fracture cases requiring mechanical traction.
 - Single patient use pad provides ease of use and reduced risk of cross contamination.
 - 23 cm x 15 cm
- Not made with natural rubber latex

Peg Board Covers

PC100

Case Pack: 6 Sets of 6 Covers



- Single patient use set of peg board covers for use with peg boards during lateral positioning.
- Soft foam with appropriate density provides protection and reduced pressure to patient.

Not made with natural rubber latex

Prime Medical, LLC's Global Headquarters is located in Knoxville, TN, and all of our products are proudly made in the USA.





CLINICALLY DRIVEN INNOVATION

Prime Medical, LLC
1410 E Weisgarber
Knoxville, TN 37909
877-434-2021

No Prime Medical, LLC products are
made with natural rubber latex.



To Order:

Email : orders@primemedicalllc.com
Fax : 877-434-2021

Additional Information:

Website:
www.primemedicalllc.com

Information requests:
info@primemedicalllc.com

Customer Service:
customerservice@primemedicalllc.com
or 877-434-2021, prompt 2