




Sterion® Vessel Loops

- Used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures
- Extruded from medical-grade silicone
- Radiopaque
- Intended for single-use and not made with natural rubber latex

Sterion® Vessel Loops

2 Loops/Pkg; 10 Pkgs/Box/Double Pouched

Part No.	Color	Cross Section – mm	Length – mm
011001PBX	Red Mini	1.5 x 0.80	406
011002PBX	Blue Mini	1.5 x 0.80	406
011003PBX	Yellow Mini	1.5 x 0.80	406
011004PBX	White Mini	1.5 x 0.80	406
011011PBX	Red Maxi	2.5 x 1.0	406
011012PBX	Blue Maxi	2.5 x 1.0	406
011013PBX	Yellow Maxi	2.5 x 1.0	406
011014PBX	White Maxi	2.5 x 1.0	406
011021PBX	White Supermaxi	5 x 1.3	559
011022PBX	Blue Supermaxi	5 x 1.3	559
011032PBX	Blue Micro	0.53 x 0.53	406
011301PBX	Red Mini	1.5 x 0.80	762
011302PBX	Blue Mini	1.5 x 0.80	762
011303PBX	Yellow Mini	1.5 x 0.80	762
011304PBX	White Mini	1.5 x 0.80	762
011311PBX	Red Maxi	2.5 x 1.0	762
011312PBX	Blue Maxi	2.5 x 1.0	762
011313PBX	Yellow Maxi	2.5 x 1.0	762
011314PBX	White Maxi	2.5 x 1.0	762