



## Technical Data Sheet Securline® Lab Marker

<b>Product Number/ Description</b>	<b>1400-16-PDC</b> Securline Red Lab Marker w/ Fine Tip <b>1400-20-PDC</b> Securline Black Lab marker w/ Fine Tip
<b>Intended Use</b>	Aspen's lab markers are versatile and can be used to make permanent markings on plastic, metal, and porcelain labware. The unique ink has the ability to write on labware even when the surface is damp.
<b>Materials of Construct</b>	The materials consist of: <ul style="list-style-type: none"><li>• Barrel, Cap and Plug- Polypropylene</li><li>• Reservoir and Nib- Polyester</li><li>• Lab Marker Ink- See attached MSDS</li></ul>
<b>Manufacturer Name and Address</b>	Aspen Surgical Products, Inc. 6945 Southbelt Dr. SE Caledonia, MI 49316 USA
<b>Legal Address of Representative of the European Union</b>	Aspen Medical Europe Ltd. Thornhill Road Redditch, Worcestershire B98 9NL England
<b>Storage</b>	No special storage requirements.
<b>Type of Sterilization</b>	N/A- Lab Markers are non-sterile.
<b>European Union Rule/ Classification / Conformity Route</b>	Lab Markers are not considered a medical device.
<b>FDA Regulatory Class and Product Code</b>	Lab Markers are not considered a medical device.